WCS®

EN Wound Care Shoe
DE Verband- und Wundtherapieschuh
ES Zapato para tratamiento de úlceras
FR Chaussure pour le soin des plaies diabétiques
IT Calzatura ortopedica terapeutica
NL Wondbehandelingsschoen
DA Bandage- og sårsko
SV Sårbehandlingssko
PL But pooperacyjny i do terapii ran
SL Čevelj za nego ran
HR Cipelja za liječenje rana
PT Sapato para terapia de feridas

DARCO

1

2


3


4


5

6


7

A


B


C


D

E

F

WCS® | Wound Care Shoe

Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This aid assists recovery and should be used exclusively under the instruction of and in consultation with your

doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product with injured, irritated, damaged or intact skin. You should use bandages or socks, for example.

We recommend that you do not drive a vehicle while you are wearing this aid. Use extreme caution when climbing stairs (e.g. also stepping up onto curbs) and walking on uneven surfaces. Walk more slowly than usual.

This aid is intended for single-patient use only.

Purpose

The Wound Care Shoe is mainly used for treating open and closed pressure sores and ulcerations (malum perforans). This enables specific offloading and pressure redistribution of the foot.

Indications / Contraindications

Indications

- > For treatment of diabetic foot syndrome
- > For healing of open wounds and and pressure ulcers and ulcerations
- > For pressure relief
- > During the postoperative phase after toe and/or forefoot surgery

- > Ulcerations up to Wagner grade 3
- > Bony prominences

Contraindications

- > Open fractures
- > Ulcerations higher than Wagner grade 3

Design features

The special sole construction of the Wound Care Shoe achieves optimum pressure redistribution. The multi-layer insole can be adapted individually to the patient to achieve selective offloading. The upper material is as seamless as possible to prevent pressing. It can be worn as an everyday shoe both indoors and outdoors.

Technical data

7 sizes:

XS (35–36), S (37–38), S/M (39–40), M (41–42), L (43–44), XL (45–46), XXL (47–48)

Universal fit for right and left

Material

Outsole: PU

Upper material: Leather

Colour: Black

Combination with other products

WCS®/WCS® Light Insoles
TwinShoe®

Maintenance

The product does not require maintenance.

Modifications may only be carried out by specialised medical personnel.

Determining the correct shoe size is very important for the functionality of the shoe.

1. Determine the correct shoe size from the insole length.
2. Open the shoe. The hook and loop fastener allows the Healing Shoe to be opened and closed quickly and easily.
3. The tried and tested 4-layer sole system is part of the delivery (Fig. 2) and can be processed individually (Fig. 6). All four insoles must always be in the shoe.
4. Position the foot with socks or bandage comfortably in the shoe.
5. A "D" is marked on the inside of the sole which should lie just behind the distal head of the metatarsal I joint. This ensures that the pressure point of the rocker-style sole is proximal to the metatarsal head.
6. Close the hook and loop straps so the foot is firmly secured.

Ulcer marking and cutting the soles:

1. Place a sterile film on the wound or pressure ulcer (Fig. 3) and circle the affected area on the film with the marker pen provided (Fig. 4).
2. The patient should now stand up on the insole in the shoe(s) so the ulcer marks the upper layer of the insole.

3. The next two layers are now punctured through the top layer using a pointed object, thus defining the cutout to be formed (Fig. 5).
4. Now draw the cutouts to be removed on the underside of the layers, as described in Fig. 7.

Leather modification (non-load-bearing sores):

1. To remove leather without damaging the liner, pass a blunt object (i.e. paperclip) from one perforated hole to the next. Be careful not to puncture the liner.
2. Cut the leather using the blunt object as a guide.
3. Insert a flat instrument between the leather and liner to continue the cut. Only cut within the given perforation holes.

Sole cut:

- A No plantar open sores
- B Severity 0 and 1, Wagner scale, open wound, cutting angle around 30°
- C After wound healing (from B)
- D Severity of 2 and higher, Wagner scale, open wound
- E After wound healing (from D)
- F Transverse section of the relief of a pressure ulcer under the third metatarsal head

Plantar layer system:

	Material	Density	Thickness	Durometer*
	Plastazote®	soft	6.3 mm	26
	EVA	medium	6.3 mm	42
	EVA	firm	6.3 mm	50
	Insole foam	high	3.1 mm	29

* The insole durometers may vary within 2° to 3° of those listed.

Note

Modifications may only be carried out by specialised medical personnel. The period of use depends on the grade of the ulcer and should be coordinated with the treating doctor. For a physiologically correct gait, use the shoe in pairs if possible or with the TwinShoe® balance shoe for height balancing.

Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling/cleaning.

Cleaning instructions

If the Healing Shoe is slightly soiled, you can wipe it with a damp cloth or clean it with a soft brush. Please do not use any aggressive cleaning agents (alcohol etc.). Do not put the shoes in the dryer or on a heat source, but leave them to dry at room temperature. The soles provided can be cleaned with lukewarm water and soap or with conventional disinfectants.

Disposal

Dispose of after use. Please dispose of the packaging in an environmentally friendly manner. The product can be disposed of with household waste.

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A. Algeo Ltd
DARCO International, Inc.