Slimflex Orthotic Insoles

Intended Purpose: Slimflex brand orthotics have been designed and developed to be simple, adaptable, vibrant and cost effective and can provide simple solutions for musculoskeletal conditions of the lower limb including overuse injuries, metatarsalgia, Morton's neuroma, heel pain, tendinopathies and osteoarthritis. They provide support, stability and pain relief and for many clinicians they are the insole of choice for a multitude of foot conditions including peripheral neuropathy and falls prevention.

Indications:

- Lower limb overuse injuries
- Acute and chronic musculoskeletal conditions of lower limb
- Hallux limitus/rigidus
- Morton's Neuroma
- Metatarsalgia
- Pes planovalgus deformity
- Adult acquired flat foot deformity
- Tibialis posterior tendon dysfunction
- Pes carovarum deformity
- Plantar heel pain
- Plantar fasciitis or plantar fasciopathy
- Achilles tendinopathy
- Generalised foot, ankle, knee or hip pain.

Contra-indications: The Slimflex range of insoles are not recommended for use in patients with diabetes unless under the explicit direction and supervision of a qualified healthcare professional. Use of the Slimflex range of inserts is not advised on areas of broken or ulcerated skin unless under the explicit direction and supervision of a qualified healthcare professional.

Serious Incident Reporting: In the event of a serious incident related to this device, report it to Algeos and the competent authority in your country (e.g. MHRA in the UK via <http://yellowcard.nhs.uk>).

Target Users: Intended for use by adults and children over the age of 5. May be modified for specific users using components specifically designed for compatibility with the Slimflex range of insoles, under the appropriate supervision and care of a qualified healthcare professional.

Lay User Guidance: If pain persists after wearing the insoles for a week, or if blisters or redness appear, discontinue use and consult a podiatrist or healthcare professional. These insoles have not been clinically tested or certified for use by individuals with diabetes. If you have altered sensation and/or circulation issues, or are at risk of developing diabetic foot ulcers, please speak with your podiatrist or diabetic foot care specialist before using them. People with diabetes may become more susceptible to altered pressure and shear forces. They may lack sufficient sensation to detect or sufficient circulation to resolve a problem which may occur as a consequence of using insoles - see residual risks section. Do not attempt to modify this device unless you are appropriately qualified. Unauthorised or improper modifications may compromise the device's safety, functionality, and compliance with regulatory standards. If you are uncertain about its operation or have any concerns, please consult a qualified healthcare professional before use.

Expected Clinical Benefits: The Slimflex brand insoles offer several clinical benefits including:

- Provides pressure redistribution to aid in pain relief, enhances comfort, and reduces fatigue during use
- Provides support to enhance stability during use.
- Offers versatility by being suitable for use in its original manufactured form or through modification with compatible supplementary components.
- Customisable to accommodate individual patient needs through modification and integration of appropriate components.
- Constructed from high-quality materials and engineered to endure routine daily use, providing reliable and long-lasting support. This product does not contain latex.
- Designed to fit a wide range of footwear types and may be transferred between different pairs as needed.
- Indicated for a broad range of applications, including everyday use and participation in sporting activities.
- Designed to enhance functional performance and assist in reducing the risk of injury during use.
- Supports improvements in gait by influencing the kinetics and biomechanics of the foot and lower limb.
- Assists in reduction and redistribution of plantar pressure, contributing to the prevention of foot ulceration.
- Contributes to enhanced quality of life by supporting improved comfort, mobility, and functional outcomes.

Performance Characteristics:

- Mouldable, semi-rigid or firm EVA base.
- Available in full length, $\frac{3}{4}$ length, and various densities.
- Colour-coded options for clinical clarity.
- Can be customised or modified for bespoke use.

Suitability & Accessory Selection:

- Select based on foot length, arch type, and clinical need.
- Accessories may include wedges, heel lifts or top covers.
- Compatible adhesives: Rubber adhesive, Polypropylene adhesive.
- Non-compatible adhesives: cyanoacrylate-based glues may damage EVA.

Residual Risks / Side Effects:

- Residual risks may include skin irritation resulting from prolonged use and discomfort due to improper fit.
- Patients should discontinue use immediately if discomfort or skin irritation occurs and consult a healthcare professional if symptoms persist.
- Side effects are uncommon; however, improper fitting may result in friction blisters or the development of pressure points.

User Specifications:

- Product colours represent density/type (e.g., dark blue = firm, light blue = soft).
- Choose appropriate size and density based on the patient profile.
- Trimmable with scissors to fit various shoe shapes.

Preparatory Handling: Insoles may be modified by a Healthcare Professional prior to final use through the addition of compatible components such as metatarsal domes, wedges, or top covers, as clinically indicated. Suitable for grinding modifications by a qualified Healthcare Professional.

Special Requirements for Use:

- Professional fitting is recommended.
 - No formal training required for general use, but professionals should be qualified in orthotic fitting and biomechanics.
 - Inappropriate or unqualified modification may impact function and increase risk of injury.
 - The insoles are not intended for continuous use however, they may be worn for periods of up to 12 hours. A gradual break-in period is recommended, with shorter durations of wear advised during the first 14 days. Always follow the guidance of a qualified healthcare professional.
-

Verification of Correct Use:

- Fit assessment: The insole should be positioned flush within the shoe, ensuring proper alignment and adequate support of the medial arch.
- Insoles should be kept clean and dry. If they become wet, allow them to dry completely before reusing to maintain performance and hygiene.
- It is recommended that insoles be worn with socks or hosiery to minimise direct skin contact and help preserve the quality of the device.
- Replace the insole if the material shows signs of deterioration or if structural integrity is compromised.
- Calibration is not required prior to use. Goods are supplied ready for immediate application, allowing them to be used directly from the packaging without any additional setup or adjustment.

Compatibility with Other Devices/Equipment:

- May be used with compatible orthotic additions (e.g., wedges, heel lifts).
- Avoid combining with adhesives or parts from other brands without confirmation of compatibility.

Warnings, Precautions, and Limitations of Use:

- Discontinue if pain worsens or if device is visibly damaged.
- Do not use if signs of excessive wear, cracking, or compression flattening appear.
- Lay users should consult a qualified healthcare professional if uncertain about fit or appropriate use of device.

Storage/Handling Conditions: Store in a cool, dry place away from direct sunlight. Avoid prolonged exposure to high heat or excessive moisture. Do not bend or distort the insole before fitting.

Care for your insole: Wipe clean with a damp cloth and mild soap. Keep away from heat. Store in a cool, dry place. This device is single patient use only and does not need to be cleaned between uses.

Disposal: Please dispose of the packaging in an environmentally friendly manner. Dispose of this product in accordance with applicable local laws and regulations.

Device Name / Trade Name: Slimflex Orthotic Insoles – available in a variety of configurations.

Registered Trademark: Slimflex.

Manufacturer Details:

A. Algeo Ltd.
Sheridan House, Bridge Industrial Estate, Speke Hall Road, Liverpool, L24 9HB (UK).
Website: www.algeos.com
Email: info@algeos.com

EU Authorised Representative (EUAR):

Advena Limited
Tower Business Centre, 2nd Floor Tower Street, Swatar, BKR 4013 (Malta).
Phone: +356 2546 6689
Email: info@advena.lt

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