

Bodytonix Elastic Knee 4-Way Support with Patella Pad

Technical characteristics and usage advice

A four way stretch knee support with silicone patella ring to help with muscle strain and knee joint pain

Features & Benefits:

- Closed kneecap
- Fits right or left leg.
- 4-way stretchable knit gives the knee comfortable compression.
- Anatomical design offers a natural fit.
- Gel patella buttress exerts pressure and stabilises the kneecap.

Technical Characteristics:

- Available in Small, Medium, Large and X-Large.
- Elasticated fabric.
- Latex free.

Indications for Use:

- Mild to moderate injuries that require support while reducing pressure to the kneecap.
- Support for weak knees, mild ligament damage and pain relief from meniscus cartilage tear.
- For use during sport activities.

Cautions:

- Its use should be supervised by a physician or health care specialist, who should provide direction on frequency and duration of use.
- In case of any abnormal sensation such as tingling, burning or numbness, pain, swelling or skin irritation, discontinue use and seek medical assistance.

Storage:

Store at room temperature. Do not expose to direct sunlight.

Markings:

Brand label.

Care Instructions:

Can be wiped clean with a damp cloth using a mild soap or disinfectant. Do not dry on a radiator.

Warranty:

12-month cover from purchase for defects in materials and workmanship.



Code	Size	Thigh Circumference
BT0648	Small	32cm
BT0650	Medium	34cm
BT0652	Large	35.5cm
BT0654	X-Large	37cm