

# Bodytonix Rigid Finger Splints

## Technical characteristics and usage advice

The Bodytonix rigid splints provide protection to the finger following injury.

### Features & Benefits:

- Finger splint protects and supports the distal joint of the finger in extension
- Ventilated design for great comfort.
- Hardwearing.
- Wipeable and reusable.
- Available in pack of 10 splints.

### Indications for Use:

- Aid in the repair of ruptured tendons and protection of the digits.

### Technical Characteristics:

- Available in 7 sizes
- Materials: high quality rigid plastic.

### Cautions:

- Its use should be supervised by a physician or health care specialist, who should provide direction on frequency and duration of use.
- In case of any abnormal sensation such as tingling, burning or numbness, pain, swelling or skin irritation, discontinue use and seek medical assistance.

### Storage:

Store at room temperature. Do not expose to direct sunlight.

### Markings:

This product contains no additional markings other than the brand label.

### Care Instructions:

Can be wiped clean with a damp cloth using a mild soap or disinfectant. Do not dry on a radiator.

### Warranty:

12-month cover from purchase for defects in materials and workmanship.



Code	Size	Length
BT0688	Size 1	4.5cm
BT0690	Size 2	5cm
BT0692	Size 3	5.5cm
BT0694	Size 4	6cm

Code	Size	Length
BT0696	Size 5	7cm
BT0698	Size 5.5	7.5cm
BT0700	Size 6	8cm
BT0702	Size 7	9cm