



SAFE WORKING LOAD

152kg (24 Stone)

PRESSURE CARE CATEGORY

High Risk

SERVICE LIFE

24-36 Months**Product Description**

Caversham® Élysée is a static dual-layer foam pressure-redistribution cushion designed for patients assessed as High Risk of pressure-ulcer development. It features a viscoelastic PU top layer for immersion and envelopment, supported by a CMHR (Combustion-Modified High Resilience) foam base to prevent bottoming out and maintain stable posture.

The waterproof, vapour-permeable, wipe-clean stretch cover (Blue) and non-slip base make Élysée suitable for wheelchairs, hospital seating, care homes and community equipment services.

Key Features

- Dual-layer construction: viscoelastic top layer + CMHR base layer
- Immersion and envelopment for High-Risk pressure care
- Stable seating geometry for posture and transfers
- Waterproof, vapour-permeable, wipe-clean stretch cover (Blue)
- Non-slip base
- Safe Working Load (SWL): 152 kg (24 stone)

Cushion Core

- Top layer: Viscoelastic PU foam
 - Density: 55–60 kg/m³
 - Firmness: ILD @40%: 100 N
- Base layer: CMHR foam
 - Density: 38–40 kg/m³
 - Firmness: ILD @40%: 200 N

Dimensions and Weight

- Size: 43 cm × 43 cm × 8.5 cm
- Cushion weight: 1.0 kg
- Safe Working Load: 152 kg (24 stone)

Cover

- Material: Waterproof, vapour-permeable, four-way stretch PU-coated textile
- Colour: Blue
- MVTR (Moisture Vapour Transmission Rate): MVTR ≥ 300 g/m²/24h.
- Hydrostatic head: 1,500 mm
- Seams Protected-seam construction
- Base: Non-slip PU underside
- **Cleaning compatibility:**
 - Compatible with standard NHS detergent wipes
 - Compatible with neutral-pH detergent cleaning
 - Compatible with chlorine-based disinfectants up to 1,000 ppm
 - Machine washable up to 60°C
 - Tumble dry on a cool setting
 - Do NOT wash above 60°C
 - Do NOT iron
 - Do NOT autoclave
 - Internal components must not be washed, soaked, submerged or machine processed

Performance & Durability

- Pressure Care Category: High - Very High Risk
- Designed for repeated wipe-clean and disinfection cycles
- Service Life: 24–36 months
- Latex-free, phthalate-free

Compliance & Regulatory

- UKCA and CE marked
- MRHA registered Class 1 Medical Device
- Fire Rating: BS 7176 – Medium Hazard
- Manufactured under ISO 9001-certified QMS
- Materials compliant with ISO 10993

Distributed by



A. Algeos Ltd
Sheridan House
Bridge Industrial Estate
Speke Hall Road
Liverpool
L24 9HB, UK

+44 (0)151 448 1228
+44 (0)151 448 1008
algeos.com
info@algeos.com