

Silverfox 5 motor Electric Podiatry Chair

with Medical Leather

SILVERFOX

Robust Electric and hydraulic chair with 5 motors and medical grade upholstery for the ultimate patient comfort.

Product Features & Benefits:

- Up & down movement with low entry height of 58 cm for wheel chair transfer and 100 cm at maximum height when seat rest tilted up to 10°. Height ranges from 58cm to 88cm when the chair is in a flat position.
- 200kg/ 440 lb safe working load tested.
- World leading LINAK 'High Reliability' actuators
- Back rest range up to 85 degrees
- The seat rest can be tilted at an angle and swirl aside to a max. of 60 degrees manually
- Leg rest area swirl aside to a max. of 90 degrees. Each leg rest swirl to a 70 degree angle.
- Seat and back controls included as standard, via hand controller.
- Cushioned plug-in headrest included as standard.
- CE Marked and Class I Medical Device, MHRA registered in the UK.
- 1 year warranty. Spare parts available.

Dimensions: L 192cm x W 80 cm (including arm rests) x H 58 cm (minimum)

Net weight: 82 kg/ 181 lb

Upholstery:

Upholstered in durable and comfortable Medical grade synthetic leather, fire retardant quality by UK Standards. Can be disinfected with alcohol or hydrogen peroxide-based daily cleaning products.

Colours: Navy blue as standard. Please contact us for other colours.



Delivery: 2-man delivery by Algeos Team, geographical location charges may apply.

Servicing/Repairs: carried out by Algeos Servicing Team.

Finance: To discuss financing options for this product, please call us on 0151 4481228.



Technical Specifications:



| Code | Description |
|--------|-----------------------------------|
| NV5304 | Silverfox Electric Podiatry Chair |



Financing Options Available



A. Algeo Ltd
 Sheridan House
 Bridge Industrial Estate
 Speke Hall Road
 Liverpool
 L24 9HB, UK

+44 (0)151 448 1228
 +44 (0)151 448 1008
 algeos.com
 info@algeos.com