

# Podotech EVA Wedge Strips

## Technical characteristics and usage advice

When it comes to making bespoke insoles, we stock a wide range of Podotech Components that provide the clinician with as many options as possible in order to attain an optimal fit.

### Features & Benefits:

- Wedged Strips made from High density (A65) EVA which have been pre-cut into wedge strips making them suitable for insole adaptations.
- Adhesive backed.
- Colour coded for easy recognition.
- Easily trimmed to be added to modular insoles.
- Sold individually.

### Technical characteristics:

**Materials:** High density (A65) EVA

**Size:** 95cm x 8cm.

**Colours:** Grey or Beige depending on degree (see table)

### Indications:

Fabrication of insole adaptations.

### Care/ Storage Instructions:

Store at room temperature. Do not expose to direct sunlight, heat or moisture.

### Markings:

This product does not have markings.

### Warranty:

12-month cover from purchase for defects in materials and workmanship.



Code	Size
OG4013	Podotech EVA Wedge Strip 3.5° Grey
OG4015	Podotech EVA Wedge Strip 5° Beige