

Podotech EVA Heel Lifts

Technical characteristics and usage advice

Features & Benefits:

- Diecut heel raises manufactured from medium-density EVA.
- A preformed adjustable wedge designed to elevate the heel.
- A linear plane, ideal for modification.
- Universal left or right.
- Sufficiently durable to resist compression.
- Fits easily into regular shoes.
- Graduated profile and self adhesive backing.
- A50 shore
- Sold in pairs

Technical characteristics:

Materials: High density (A50) EVA

Size: available in 3 sizes

Colours: Blue

Indications:

Fabrication of insole adaptations.

Care/ Storage Instructions:

Store at room temperature. Do not expose to direct sunlight, heat or moisture.

Markings:

This product does not have markings.

Warranty:

12-month cover from purchase for defects in materials and workmanship.



Heel Lift Size	Dimensions	Typical UK size
Small	L 7.5cm x W 5.5cm	2 to 5
Medium	L 8.5cm x W 6cm	6 to 8
Large	L 9.5cm x W 6.8cm	9 to 14