







A. Algeo Ltd
 Sheridan House
 Bridge Industrial Estate
 Speke Hall Road
 Liverpool
 L24 9HB, UK








+44 (0)151 448 1228
 +44 (0)151 448 1008
 algeos.com
 info@algeos.com








Product Size	Shoe UK Size Recommended
Small (s. 1)	2 to 4.5
Medium (s.2)	5 to 9.5
Large (s.3)	10 and above





Metatarsal and Valgus Pads

Note dimensions are subject to batch variations

Product Image	Order Code	Description	Length (mm)	Width (mm)	Height (mm)
	SD2001	Metatarsal Pads - Medium Firm Rubber - Size: 1 x 10 Pairs	48	38	3
	SD2002	Metatarsal Pads - Medium Firm Rubber - Size: 2 x 10 Pairs	51	42	4
	SD2003	Metatarsal Pads - Medium Firm Rubber - Size: 3 x 10 Pairs	56	39	5
	SD2021	Metatarsal Pads - High Soft Rubber - Size: 1 x 10 Pairs	51	36	3
	SD2022	Metatarsal Pads - High Soft Rubber - Size: 2 x 10 Pairs	55	38	3
	SD2023	Metatarsal Pads - High Soft Rubber - Size: 3 x 10 Pairs	58	40	3
	SD2071	Metatarsal Pads - High Soft Foam - Size: 1 x 10 Pairs	50	32	3
	SD2072	Metatarsal Pads - High Soft Foam - Size: 2 x 10 Pairs	55	36	5
	SD2073	Metatarsal Pads - High Soft Foam - Size: 3 x 10 Pairs	62	39	8
	SD2080	Metatarsal Pads - Flat Soft Foam - Size: 1 x 10 Pairs	64	42	2
	SD2081	Metatarsal Pads - Flat Soft Foam - Size: 2 x 10 Pairs	67	45	2
	SD2091	Child's Metatarsal Pads - Size: 1 x 10 Pairs	44	32	3
	SD2092	Child's Metatarsal Pads - Size: 2 x 10 Pairs	50	35	3.5
	SD2093	Child's Metatarsal Pads - Size: 3 x 10 Pairs	55	40	4

	SD2121	Metatarsal Pads - High & Wide - Size: 1/ Small x 10 Pairs	46	38	2
	SD2122	Metatarsal Pads - High & Wide - Size: 2/ Small x 10 Pairs	50	40	3
	SD2123	Metatarsal Pads - High & Wide - Size: 3/ Medium x 10 Pairs	54	45	4
	SD2124	Metatarsal Pads - High & Wide - Size: 4/Medium x 10 Pairs	56	50	5
	SD2125	Metatarsal Pads - High & Wide - Size: 5/Large x 10 Pairs	60	50	6
	SD2126	Metatarsal Pads - High & Wide - Size: 6/Large x 10 Pairs	62	52	7
	SD2271	Metatarsal Pads - Medical Blue Foam - Small x 10 Pairs	50	45	4
	SD2272	Metatarsal Pads - Medical Blue Foam - Medium x 10 Pairs	54	45	6
	SD2273	Metatarsal Pads - Medical Blue Foam - Large x 10 Pairs	64	55	7
	SD3033	SBR Meta Pad Kidney Shape with Central Fitting - Size: 1 x 10 Pairs	73	51	5
	SD3034	SBR Meta Pad Kidney Shape with Central Fitting - Size: 2 x 10 Pairs	74	52	5
	SD3035	SBR Meta Pad Kidney Shape with Central Fitting - Size: 3 x 10 Pairs	75	53	5
	SD3036	SBR Meta Pad - Drop Shape with High Central Lift - One Size x 10 Pairs	50	40	1
	SD3039	SBR Meta Pad - Stomach Shape - Size: 1 x 10 Pairs	70	50	6
	SD3040	SBR Meta Pad - Stomach Shape - Size: 2 x 10 Pairs	71	52	6
	SD3041	SBR Meta Pad - Wide Fit - One Size x 10 Pairs	55	50	5
	SD3042	SBR Meta Pad - Kidney Shape - Size: 1 x 10 Pairs	72	35	4
	SD3043	SBR Meta Pad - Kidney Shape - Size: 2 x 10 Pairs	75	37	4.5
	SD3044	SBR Meta Pad - Kidney Shape - Size: 3 x 10 Pairs	80	39	5
	SD3045	SBR Meta Pad - Kidney Shape - Size: 4 x 10 Pairs	87	42	5

	SD3050	SBR Toe Bar - One Size x 10 Pairs	62	23	2
	SD3051	SBR Meta Pad - Drop Shape - Size: 1 x 10 Pairs	50	32	3
	SD3052	SBR Meta Pad - Drop Shape - Size: 3 x 10 Pairs	51	39	3
	SD3066	SBR Central Meta Pad - Size: 1 x 10 Pairs	88	37	4
	SD3067	SBR Central Meta Pad - Size: 2 x 10 Pairs	100	44	5.5
	SD3068	SBR Central Meta Pad - Size: 3 x 10 Pairs	110	50	6
	SD3071	SBR Cobra Pad - High Internal Arch - UK 3-5 x 5 Pairs	165	70	1
	SD3072	SBR Cobra Pad - High Internal Arch - UK 6-7 x 5 Pairs	180	70	1
	SD3073	SBR Cobra Pad - High Internal Arch - UK 8-9 x 5 Pairs	190	72	1
	SD3074	SBR Cobra Pad - High Internal Arch - UK 10+ x 5 Pairs	205	80	1
	SD3075	Valgus Wedge Pad - Thickness: 5mm - Size: 1 x 5 Pairs	70	35	3
	SD3076	Valgus Wedge Pad - Thickness: 5mm - Size: 2 x 5 Pairs	80	45	3
	SD3079	Wide Over - Arch Pad - Thickness: 5mm - One Size x 5 Pairs	92	44	3
	SD3080	Sponge Rubber Meta Pad - Drop Shape - Size: 1 x 10 Pairs	59	48	2.5

	SB1130	Soft Valgus Pads - Grey - Size: 1/Small x 5 Pairs	55	30	6
	SB1132	Soft Valgus Pads - Grey - Size: 2/Small x 5 Pairs	68	32	6
	SB1134	Soft Valgus Pads - Grey - Size: 3/Medium x 5 Pairs	85	40	6
	SB1136	Soft Valgus Pads - Grey - Size: 4/Medium x 5 Pairs	78	52	7
	SB1138	Soft Valgus Pads - Grey - Size: 5/Large x 5 Pairs	100	60	8.5
	SB1140	Soft Valgus Pads - Grey - Size: 6/Large x 5 Pairs	128	66	9
	SB1050	Valgus Pads - Med Firm, High Raise - Size: 0/Small x 5 Pairs	80	40	2
	SB1051	Valgus Pads - Med Firm, High Raise - Size: 1/Small x 5 Pairs	85	45	2
	SB1052	Valgus Pads - Med Firm, High Raise - Size: 2/Medium x 5 Pairs	95	50	2.5
	SB1053	Valgus Pads - Med Firm, High Raise - Size: 3/Large x 5 Pairs	105	50	3
	SB1102	Valgus Pads - Black - Medium x 5 Pairs	80	45	8
	SB1103	Valgus Pads - Black - Large x 5 Pairs	92	52	9
	SB1104	Valgus Pads - Black - X Large x 5 Pairs	120	57	12
	SE1181	Poron Cobra Pads - Size: 1 UK 3-4 x 5 Pairs	15	58	3
	SE1182	Poron Cobra Pads - Size: 2 UK 5-6 x 5 Pairs	16	58	3
	SE1183	Poron Cobra Pads - Size: 3 UK 8-9 x 5 Pairs	17	59	3
	SE1184	Poron Cobra Pads - Size: 4 UK 10 x 5 Pairs	18	60	3
	SE1185	Poron Cobra Pads - Size: 5 UK 10+ x 5 Pairs	19.5	60	3