

The Algeos 3D top-cover range has been developed to give clinicians flexible, activity-specific options that enhance comfort, performance and durability. Each cover is made from high-quality synthetic microfibre with integrated Poron[®] cushioning, temperature-regulating properties and excellent moisture management.



Designed for use with the XFix fastening system, all covers attach securely without glue and can be replaced or upgraded in seconds, supporting cleaner workflows and long-term patient adaptability.

Key Features

Interchangeable with XFix:

All top covers attach without adhesives, allowing quick replacement and reducing waste.

High-Performance Microfibre Construction:

Durable, breathable and soft against the skin, suitable for daily and active use.

Poron[®] Cushioning:

Delivers long-lasting shock absorption, comfort and resilience.

Temperature-Regulating, Antibacterial & Antimicrobial:

Helps maintain dryness and hygiene across long periods of wear.

Washable & Hypoallergenic:

Helps maintain dryness and hygiene across long periods of wear.

Range-Specific Performance:

Includes options for sport propulsion, daily comfort, impact protection and targeted clinical offloading.



Everyday Specification

Name	Everyday: Leathertec + PORON® Performance 4000
Texture	Synthetic microfibre
Colour	Charcoal
Material	100% microfibre
Base	Velour fastening
Thickness	3.2 mm
Density	20 ft ³
Breathability	Absorption–desorption = complete dryness
Shock absorption	Medium
Patient comfort	Temperature-regulating / antibacterial / antimicrobial
Touch	Smooth
Sweat absorption	Sweat-resistant
Hypoallergenic	Yes
Maintenance	Washable (keeps all properties)
Recommended for	Trainers / Business / Casual / Outdoor
Brandsole options	Standard
Size range	35–55



Sport Specification

Name	Sport: Leathertec + PORON® Vive
Texture	Synthetic microfibre
Colour	Charcoal
Material	100% microfibre
Base	Velour fastening
Thickness	3.2 mm
Density	16 ft ³
Breathability	Complete dryness
Shock absorption	High propulsion
Patient comfort	Temperature-regulating/antibacterial/antimicrobial
Patient comfort	Smooth
Sweat absorption	Sweat-resistant
Hypoallergenic	Yes
Maintenance	Washable (keeps all properties)
Recommended for	Trainers
Brandsole options	Standard
Size range	35–55



Protect Specification

Name	
Texture	
Colour	Charcoal
Material	100% microfibre
Base	Velour fastening
Thickness	3.2 mm
Density	15 ft ³
Breathability	Absorption–desorption = complete dryness
Shock absorption	High
Patient comfort	Temperature-regulating / antibacterial / antimicrobial
Touch	Smooth
Sweat absorption	Sweat-resistant
Hypoallergenic	Yes
Maintenance	Washable (keeps all properties)
Recommended for	Trainers / Business / Casual / Outdoor
Brandsole options	Standard
Size range	35–55



Care Hindfoot Specification

Name	Care Hindfoot: Leathertec + PORON® Medical in Hindfoot Zone
Texture	Synthetic microfibre
Colour	Charcoal
Material	100% microfibre
Base	Velour fastening
Thickness	3.2 mm
Density	15 ft ³ heel / 20 ft ³ forefoot
Breathability	Absorption–desorption = complete dryness
Shock absorption	Low–Medium
Patient comfort	Temperature-regulating / antibacterial / antimicrobial
Touch	Smooth
Sweat absorption	Sweat-resistant
Hypoallergenic	Yes
Maintenance	Washable (keeps all properties)
Recommended for	Business / Casual
Brandsole options	Standard
Size range	35–55



Care Forefoot Specification

Name	
Texture	
Colour	Charcoal
Material	100% microfibre
Base	Velour fastening
Thickness	2.4 mm
Density	20 ft ³ heel / 15 ft ³ forefoot
Breathability	Absorption–desorption = complete dryness
Shock absorption	Medium–Low
Patient comfort	Temperature-regulating / antibacterial / antimicrobial
Touch	Smooth
Sweat absorption	Sweat-resistant
Hypoallergenic	Yes
Maintenance	Washable (keeps all properties)
Recommended for	Business / Casual
Brandsole options	Standard
Size range	35–55